

# Scissors, Lister Bandage

These bandage scissors feature an angled design and blunt tip for safe, controlled cutting. They are ideal for cutting bandages, dressings, and other materials during post-mortem procedures.



## SPECIFICATIONS

### Construction

- Tips: Blunt/probe-style lower tip
- Cutting Edge: Angled
- Available in Premium or Standard grades

### Materials

- Stainless steel

### Dimensions

- AA060 & AA160: 3.5" L
- AA061 & AA161: 4.5" L
- AA062 & AA162: 5.5" L
- AA063: 7.25" L

## Variations

| Model No. | Description              | Grade    | Dimensions |
|-----------|--------------------------|----------|------------|
| AA060     | Scissors, Lister Bandage | Standard | 3.5" L     |
| AA061     | Scissors, Lister Bandage | Standard | 4.5" L     |
| AA062     | Scissors, Lister Bandage | Standard | 5.5" L     |
| AA063     | Scissors, Lister Bandage | Standard | 7.25" L    |
| AA160     | Scissors, Lister Bandage | Premium  | 3.5" L     |
| AA161     | Scissors, Lister Bandage | Premium  | 4.5" L     |
| AA162     | Scissors, Lister Bandage | Premium  | 5.5" L     |

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI