


Scissors, Potts-Smith, Angled

These scissors feature an angled tip for improved access during post-mortem dissection. Their design allows for controlled cutting in hard-to-reach areas.



SPECIFICATIONS

Construction

- Tips: Sharp

Materials

- Highest quality German stainless steel

Dimensions

- AA144 & AA145: 7.5" L
- AA144: 60 degree angle
- AA145: 45 degree angle

Variations

Model No.	Description	Angle Degree	Dimensions
AA144	Scissors, Potts-Smith	60	7.5" L
AA145	Scissors, Potts-Smith	45	7.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI