

# Scissors, Iris Dissection

Also referred to as the “Ribbon Scissor” because of the unique and finger-friendly, comfort handle. With extra wide finger rings offers comfortable post-mortem dissection. This precise and delicate scissor is made of the highest quality stainless steel. Extra sharp tips allow for intricate work. Curved or straight.



## SPECIFICATIONS

### Construction

- Tips: Extra Sharp
- Available in Straight or Curved tips

### Materials

- High grade stainless steel

### Dimensions

- AA146 & AA147: 4" L

## Variations

| Model No. | Description               | Tip      | Dimensions |
|-----------|---------------------------|----------|------------|
| AA146     | Scissors, Iris Dissection | Straight | 4" L       |
| AA147     | Scissors, Iris Dissection | Curved   | 4" L       |

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI