

SCISSORS

Scissors, Delicate, Double Curved, Sharp

Finger-friendly ribbon handle. Extra sharp with single serrated edge. Double curvature of these scissors makes post-mortem dissection easier in hard to reach areas or where a clear view is not always available.



AA154

SPECIFICATIONS

Construction

- Tips: Sharp, one side of scissor edge is serrated
- Cutting Edge: Curved

Materials

- Premium grade stainless steel

Dimensions

- 4.625" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI