

Scissors, Tungsten Carbide, Iris

The unique tungsten carbide cutting edge provides effortless, yet accurate cuts through the toughest of tissue. Scissor shanks are tempered to the ideal hardness to allow proper flexibility and comfort. Finger rings are well rounded and polished allowing better feel while preventing finger fatigue during complex post-mortem dissection cases.



SPECIFICATIONS

Construction

- Tips: Sharp/sharp points
- Available in Curved or Straight cutting edge

Materials

- Premium grade stainless steel

Dimensions

- AA405 & AA406: 4.5" L

Variations

Model No.	Description	Cutting Edge
AA405	Scissors, Tungsten Carbide, Iris	Straight
AA406	Scissors, Tungsten Carbide, Iris	Curved

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI