

Scissors, Tungsten Carbide, Metzenbaum

The unique tungsten carbide cutting edge provides effortless, yet accurate cuts through the toughest of tissue. Scissor shanks are tempered to the ideal hardness to allow proper flexibility and comfort. Finger rings are well rounded and polished allowing better feel while preventing finger fatigue during complex post-mortem dissection cases.



SPECIFICATIONS

Construction

- Tips: Narrow, blunt tips
- Available in Curved or Straight cutting edge

Materials

- Tungsten Carbide

Dimensions

- AA407 & AA408: 5.5" L
- AA409 & AA410: 7" L
- AA411 & AA412: 8" L
- AA413 & AA414: 9" L
- AA415: 11" L

Variations

Model No.	Description	Cutting Edge	Dimensions
AA407	Scissors, Tungsten Carbide	Straight	5.5" L
AA408	Scissors, Tungsten Carbide	Curved	5.5" L
AA409	Scissors, Tungsten Carbide	Straight	7" L
AA410	Scissors, Tungsten Carbide	Curved	7" L
AA411	Scissors, Tungsten Carbide	Straight	8" L
AA412	Scissors, Tungsten Carbide	Curved	8" L
AA413	Scissors, Tungsten Carbide	Straight	9" L
AA414	Scissors, Tungsten Carbide	Curved	9" L
AA415	Scissors, Tungsten Carbide	Curved	11" L

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI