

Scissors, Tungsten Carbide, Mayo

The unique tungsten carbide cutting edge provides effortless, yet accurate cuts through the toughest of tissue. Scissor shanks are tempered to the ideal hardness to allow proper flexibility and comfort. Finger rings are well rounded and polished allowing better feel while preventing finger fatigue during complex post-mortem dissection cases.



SPECIFICATIONS

Construction

- Tips: Blunt rounded tips
- Available in Curved or Straight cutting edge

Materials

- Tungsten Carbide

Dimensions

- AA416 & AA417: 5.5" L
- AA418 & AA419: 6.75" L
- AA420 & AA421: 9" L

Variations

Model No.	Description	Cutting Edge	Dimensions
AA416	Scissors, Tungsten Carbide	Straight	5.5" L
AA417	Scissors, Tungsten Carbide	Curved	5.5" L
AA418	Scissors, Tungsten Carbide	Straight	6.75" L
AA419	Scissors, Tungsten Carbide	Curved	6.75" L
AA420	Scissors, Tungsten Carbide	Straight	9" L
AA421	Scissors, Tungsten Carbide	Curved	9" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI