

SCISSORS

Scissors, SuperCut, Sharp/Blunt

SuperCut scissors are tempered and hardened for an extended, sharper life over ordinary steel. The SuperCut blade has a specially designed sharpened edge. They also feature rounded and polished finger rings. A general tool for cutting.



SPECIFICATIONS

Construction

- Tips: Sharp/blunt
- Available in Curved or Straight cutting edge

Materials

- Stainless steel

Dimensions

- AA508 & AA511: 5.5" L

Variations

Model No.	Description	Cutting Edge	Dimensions
AA508	Scissors, SuperCut	Straight	5.5" L
AA511	Scissors, SuperCut	Curved	5.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI