


Scissors, SuperCut Mayo

SuperCut scissors are tempered and for a longer, sharper life over ordinary steel. The SuperCut blade has a specially designed razor sharp edge. The straight version can cut into surface area tissue while the curved version can cut into heavier or deeper tissue. They also feature rounded and polished finger rings.



SPECIFICATIONS

Construction

- Tips: Slightly rounded tips
- Available in Curved or Straight cutting edge

Materials

- Stainless steel

Dimensions

- AA519 & AA521: 5.5" L
- AA520 & AA522: 6.75" L
- AA537 & AA538: 9" L

Variations

Model No.	Description	Cutting Edge	Dimensions
AA519	Scissors, SuperCut Mayo	Straight	5.5" L
AA520	Scissors, SuperCut Mayo	Straight	6.75" L
AA521	Scissors, SuperCut Mayo	Curved	5.5" L
AA522	Scissors, SuperCut Mayo	Curved	6.75" L
AA537	Scissors, SuperCut Mayo	Straight	9" L
AA538	Scissors, SuperCut Mayo	Curved	9" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI