

Scissors, SuperCut Metzenbaum, Delicate

SuperCut scissors are tempered and hardened for a longer, sharper life over ordinary steel. The SuperCut blade has a specially designed keen edge combined with a dull tip, which can be useful for delicate post-mortem dissection. Also feature rounded and polished finger rings.



SPECIFICATIONS

Construction

- Tips: Narrow blunted tips
- Cutting Edge: Razor sharp cutting edge
- Available in Curved or Straight cutting edge

Materials

- Premium grade

Dimensions

- AA523 & AA525: 5" L

Variations

Model No.	Description	Cutting Edge	Dimensions
AA523	Scissors, Delicate	Straight	5" L
AA525	Scissors, Delicate	Curved	5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI