

Scissors, SuperCut Metzenbaum, Curved

SuperCut scissors are tempered and hardened for a longer sharper life than ordinary steel. The SuperCut blade has a specially designed sharpened edge. The scissors also feature rounded and polished finger rings. Manufactured with a curve to attain the hard to reach delicate tissue or angle post-mortem dissections.



SPECIFICATIONS

Construction

- Tips: Narrow blunt tips
- Cutting edge: Razor sharp cutting edge

Materials

- Stainless steel

Dimensions

- AA526: 7" L
- AA539: 9" L

Variations

Model No.	Description	Dimensions
AA526	Scissors, SuperCut Metzenbaum	7" L
AA539	Scissors, SuperCut Metzenbaum	9" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI