

# Scissors, SuperCut Stevens Tenotomy

SuperCut scissors are tempered and hardened for a longer, sharper life than ordinary steel. The SuperCut blade features a specially designed acute edge, while rounded and polished finger rings provide comfortable handling for delicate post-mortem dissection.



## SPECIFICATIONS

### Construction

- Tips: Narrow sharp tips
- Cutting Edge: Razor sharp cutting edge
- Available in Curved or Straight cutting edge

### Materials

- Premium grade stainless steel

### Dimensions

- AA535 & AA536: 4.5" L

## Variations

| Model No. | Description        | Cutting Edge | Dimensions |
|-----------|--------------------|--------------|------------|
| AA535     | Scissors, SuperCut | Straight     | 4.5" L     |
| AA536     | Scissors, SuperCut | Curved       | 4.5" L     |

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI