

Disposable Scissors

These disposable scissors are designed for convenient use during post-mortem procedures. They are highly recommended for infectious cases where single-use instruments are preferred.



SPECIFICATIONS

Construction

- Cutting Edge: Straight
- Available in Sharp/Sharp or Sharp/Blunt tip

Materials

- Nickel-plated carbon steel

Dimensions

- AA601 & AA602: 4.5" L
- AA604: 5.25"

Quantity

- 12 per pack

Performance

- Not autoclavable

Variations

Model No.	Description	Tip	Dimensions
AA601	Scissors, Disposable	Sharp/Sharp	4.5" L
AA602	Scissors, Disposable	Sharp/Blunt	4.5" L
AA604	Scissors, Disposable	Sharp/Blunt	5.25" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI