

Scissors, Titanium, Mayo

These Titanium Alloy Mayo scissors are designed for durable, precise cutting during post-mortem dissection. The titanium construction helps maintain sharp cutting edges longer and provides greater stain resistance than traditional stainless steel instruments.



SPECIFICATIONS

Construction

- Available in Curved or Straight cutting edge

Dimensions

- AA701 & AA702: 6.75" L

Performance

- NOTE: Do not use bleach or other caustic chemicals

Materials

- Standard grade titanium alloy

Variations

Model No.	Description	Cutting Edge	Dimensions
AA701	Scissors, Titanium, Mayo	Straight	6.75" L
AA702	Scissors, Titanium, Mayo	Curved	6.75" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI