


Scissors, Titanium, Metzenbaum

These Titanium Alloy scissors are designed for durable, precise cutting during post-mortem dissection. The titanium construction helps maintain sharp cutting edges longer and provides greater stain resistance than traditional stainless steel instruments.



SPECIFICATIONS

Construction

- Available in Curved or Straight cutting edge

Dimensions

- AA703 & AA704: 5.5" L
- AA705 & AA706: 7" L

Performance

- NOTE: Do not use bleach or other caustic chemicals

Materials

- Standard grade titanium alloy

Variations

Model No.	Description	Cutting Edge	Dimensions
AA703	Scissors, Titanium	Straight	5.5" L
AA704	Scissors, Titanium	Curved	5.5" L
AA705	Scissors, Titanium	Straight	7" L
AA706	Scissors, Titanium	Curved	7" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI