

# Scissors, Enterotomy

These sturdy, large gastrointestinal scissors are designed for controlled cutting during post-mortem dissection. Their specialized design allows tissue to be gently pulled away while cutting.



## SPECIFICATIONS

### Construction

- Tips: Hook blade
- Available in Standard or Premium grades

### Materials

- Stainless steel

### Dimensions

- AH003 & AH103: 8.25" L

## Variations

| Model No. | Description          | Grade    | Dimensions |
|-----------|----------------------|----------|------------|
| AH003     | Scissors, Enterotomy | Standard | 8.25" L    |
| AH103     | Scissors, Enterotomy | Premium  | 8.25" L    |

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI