

Spray Rinse, Control Valve and Vacuum Breaker

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



OA-215

SPECIFICATIONS

Construction

- Dual Check Vacuum Breaker with Intermediate Vent
- Control Valve with Swiveling Seat Disk
- Glass Filler Valve
- Insulated Handle

Materials

- PVC Clearbraid Tubing
- Spiraled Polyester Yarn

Dimensions

- 5/16" inner diameter and 17/32" outer diameter

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI