

End Rinse Assembly

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



OA-300

SPECIFICATIONS

Construction

- Perimeter rinse attaches to pedestal autopsy tables (not wall-mounted).
- Sits under the workstation and drains into a communal sink
- Attaches to a cold-water line
- (7) nipple spouts along the perimeter rinse

Materials

- 16-gauge 304 Stainless steel square tubing with a 180-grit finish
- HDPE nipple spouts

Dimensions

- 1.2" L X 20" W X 1.8"H

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI