



Three-Sided Rinse Assembly for Autopsy Tables

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



SPECIFICATIONS

Construction

- Perimeter rinse attaches to pedestal autopsy tables. It sits under the workstation and drains into a communal sink
- Attaches to a cold-water line
- There are small perforations along the perimeter rinse to wash any material into the sink

Materials

- 304 Stainless steel tubing with a 180-grit finish

Dimensions

- 74" L X 27" W X 1" H

Note: Photos shown are not exact product. Photos are for reference only.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI