




REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM

Item #	Qty	File Name	Title	Material	Rev.
1	2	PB0923.par	SCREW, SOCKET HEAD CAP- #8-32 X 1/4" L, SS		A
2	1	PM0956.par	RULER, 78", IN/CM, ALUM		

MOPEC LEGAL NOTE
THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER

MOPEC 2026

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES
TOLERANCES

0.0	± 0.1 in
0.00	± 0.05 in
0.000	± 0.010 in
0.0000	± 0.0010 in
ANGLES	± 0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MATERIAL

MOPEC GROUP
ADVANCED TECHNOLOGY MANUFACTURING

DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: bdevins
DATE: 2/9/2026
CHECKED BY: ---
DATE: 2/9/2026
APPROVED BY: MMcCLAIN
DATE: 2/9/2026

PART NUMBER: OA-302-A
DESCRIPTION: INCH/CM RULER, REMOVABLE, 78"
SIZE: A2
SCALE: 3:1

FILE NAME: OA-302-A.dft
REV: A

DO NOT SCALE IF IN DOUBT ASK

SHEET 1 OF 1

OA-302-A.dft