


Stainless Steel Grid Plates

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



SPECIFICATIONS

Construction

- Perforated stainless steel grid plates for pedestal autopsy tables with 0.5" perforations to allow for drainage

Materials

- 14-gauge 304 perforated stainless steel plates with a No.4 finish

Capacity

- 100 lbs per grid plate, 300 lbs total

Dimensions

- 50.8" L X 20.6" W X 1.3" D
- Each panel is 16.9" L, there are (3) panels in total
- 0.5" Perforations in a grid pattern

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI