


Adjustable Tie Down Assembly

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



SPECIFICATIONS

Construction

- Length of the tie down assembly is dependent on the length of necropsy table. Starting 5" in every 10" there is a .762 diameter whole drilled through the outer assembly. This continues until the last hole is 5" away from the end of the outer assembly. The outer assembly then screwed into the inner assembly along the drilled holes and spacers. The spacers provide a 1" space along the interior of the assembly

Materials

- 304 stainless steel assembly with a No. 4 finish

Dimensions

- 0.25" W and 1.5" D

Options & Accessories

- OA-503: Nylon Straps (used together with OA-502)

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI