


Inch/Cm Removable Ruler

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



SPECIFICATIONS

Construction

- Read in 1/16" increments on top and 1mm increments on bottom
- 48" or 72" options available

Materials

- 304 stainless steels with a No. 4 Finish
- Adhesive backed rule applied to stainless steel strip

Dimensions

- OA-504: 48" L X 1" W X 0.01" D
- OA-303: 72" L X 1" W X 0.01" D

Variations

Model No.	Description	Dimensions
OA-504	48" Inch/Cm Removable Ruler	48" L X 1" W X 0.01" D
OA-302	72" Inch/Cm Removable Ruler	72" L X 1" W X 0.01" D

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI