

Spray Hose Assembly With Backflow & Recoil Hose

Mopec's Options & Accessories are designed to enhance the functionality, efficiency, and customization of Mopec Equipment. These components provide users with flexible solutions to meet specific workflow needs, improve safety, and support specialized procedures.



OA-506

SPECIFICATIONS

Construction

- 9' Coiled polyurethane coiled hose allows for the hose to reach across most work surfaces
- Button or handle trigger allows for control over water flow

Materials

- Chrome plated brass handle
- 9' Polyurethane hose
- Brass swivel nut connection

Utilities / Power

- 5.05 GPM
- Pressure: Min: 20 psi, Max: 125 psi
- Temperature: Min: 40°F (4.4°C), Max: 140°F (60°C)

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI