

Scissors, Metzenbaum, Curved

These scissors feature slim, blunt, curved tips and thin, delicate shafts for precise post-mortem dissection with minimal damage to surrounding tissues. Available in a variety of lengths and grades to accommodate different dissection needs.



SPECIFICATIONS

Construction

- Tips: Narrow blunt points
- Cutting Edge: Curved
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA025 & AA125: 5.5" L
- AA026 & AA126: 7" L
- AA039 & AA139: 9" L
- AA040 & AA140: 11" L

Variations

Model No.	Description	Grade	Dimensions
AA025	Scissors, Metzenbaum	Standard	5.5" L
AA026	Scissors, Metzenbaum	Standard	7" L
AA039	Scissors, Metzenbaum	Standard	9" L
AA040	Scissors, Metzenbaum	Standard	11" L
AA125	Scissors, Metzenbaum	Premium	5.5" L
AA126	Scissors, Metzenbaum	Premium	7" L
AA139	Scissors, Metzenbaum	Premium	9" L
AA140	Scissors, Metzenbaum	Premium	11" L

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI