

Scissors, Metzenbaum, Straight

These scissors feature slim, blunt tips and thin, delicate shafts for precise post-mortem dissection. Their specialized design helps minimize damage to surrounding tissues while providing controlled cutting. Variety of lengths and grades.



SPECIFICATIONS

Construction

- Tips: Narrow blunt points
- Cutting Edge: Straight
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA023 & AA123: 5.5" L
- AA024 & AA124: 7" L
- AA037 & AA137: 9" L
- AA038 & AA138: 11" L

Variations

Model No.	Description	Grade	Dimensions
AA023	Scissors, Metzenbaum	Standard	5.5" L
AA024	Scissors, Metzenbaum	Standard	7" L
AA037	Scissors, Metzenbaum	Standard	9" L
AA038	Scissors, Metzenbaum	Standard	11" L
AA123	Scissors, Metzenbaum	Premium	5.5" L
AA124	Scissors, Metzenbaum	Premium	7" L
AA137	Scissors, Metzenbaum	Premium	9" L
AA138	Scissors, Metzenbaum	Premium	11" L

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI