

Transfer Board

The T3627 Transfer Board is built from high density polyethylene with 8 hand slots around the perimeter, providing a stable transfer for cadavers from cooler to table. Ideal to be used between multi-level work surfaces such as storage racks and mortuary refrigeration.



SPECIFICATIONS

Materials

- 1/2" HPDE(High density polyethylene)

Dimensions

- 23"W x 76"L x .5"H
- Approximate weight: 30lbs

Capacity

- 300lbs

Performance

- (8) hand slots with rounded edges
- Washable

Compatibility

- Mopec IJ Racks (Full length roller racks)
- Full length roller lifts
- Full length roller carts

Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
FD010	23" Transfer Board	Replaced	
FD011	30" Transfer Board	Obsolete	
FD012	36" Transfer Board	Obsolete	

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI