


Forceps, Adson Tissue

These Adson tissue forceps are designed for precise tissue handling during post-mortem dissection. The straight 1x2 toothed tips and extra-broad ribbed grip provide a secure, controlled grasp for delicate dissection procedures.



AB028

SPECIFICATIONS

Construction

- Tip: Straight, 1x2 teeth tips
- Grip: Extra broad ribbed grip

Materials

- Standard grade stainless steel

Dimensions

- 4.75" (121 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI