

SPECIALTY

Forceps, Iris Tissue, Straight, Serrated

These straight iris tissue forceps are designed for delicate tissue handling during post-mortem dissection. The serrated tips provide a secure, controlled grasp, while the ribbed grip supports precise handling during detailed dissection procedures.



AB035

SPECIFICATIONS

Construction

- Tips: Straight, serrated
- Grip: Ribbed

Materials

- Standard grade stainless steel

Dimensions

- 4" (102 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI