


SPECIALTY

Forceps, Iris Tissue, Curved

These curved iris tissue forceps are designed for precise tissue handling during post-mortem dissection. The fully curved, serrated tips provide a secure grip with minimal tissue compression, while the ribbed grip supports controlled handling.



AB036

SPECIFICATIONS

Construction

- Tips: Fully curved, serrated
- Grip: Ribbed

Materials

- Standard grade stainless steel

Dimensions

- 4" (102 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI