

Forceps, Angle Tip, Serrated

These angled-tip forceps are designed for secure, precise tissue handling during post-mortem dissection. The extra-fine serrated tips provide a firm, controlled grip, while the angled design allows improved access and manipulation in confined areas.



AB040

SPECIFICATIONS

Construction

- Tips: Extra-fine, angled, serrated
- Grip: Ribbed

Materials

- Premium grade stainless steel

Dimensions

- Overall length: 5.5"
- Handle length: 2.5"
- Forceps length from handle: 3"

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI