

SPECIALTY

Forceps, Iris Tissue, Fine Point

These fine-point iris tissue forceps are designed for precise tissue handling during post-mortem dissection. The curved, fine-point tips provide a secure, controlled grasp with minimal tissue compression, while the ribbed grip supports precise handling.



AB041

SPECIFICATIONS

Construction

- Tips: Curved, fine point
- Grip: Ribbed

Materials

- Standard grade stainless steel

Dimensions

- 4.5" (114 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI