


Forceps, Dressing, 6"

These dressing forceps are designed for handling tissue and dressing materials during post-mortem examinations and dissection procedures. The smooth tips provide controlled handling, while the ribbed grip supports a secure, precise grasp.



AB044

SPECIFICATIONS

Construction

- Tips: Smooth
- Grip: Ribbed

Materials

- Standard grade stainless steel

Dimensions

- 6" (152 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI