


Forceps, Needle Core Biopsy

These needle core biopsy forceps are designed for secure handling of core biopsy tissue during post-mortem procedures. The modified hooked tip helps capture and hold tissue samples while supporting precise, controlled handling.



AB049

SPECIFICATIONS

Construction

- Tips: Modified hooked tip
- Grip: Smooth

Materials

- Standard grade stainless steel

Dimensions

- 4.25" (108 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI