


# Rubber Tipped Forceps

These rubber-tipped forceps are designed for secure tissue handling during post-mortem dissection with minimal tissue damage or compression. The soft, rounded rubber tips prevent metal exposure and provide strong contrast for specimen photography.



## SPECIFICATIONS

### Construction

- Tips: Soft, rounded rubber
- Blue tips for specimen contrast

### Materials

- Standard grade stainless steel
- Rubber-tipped ends

### Dimensions

- Available in 6", 8", 10", or 12" L

### Performance

- Rinse under water or wipe clean with an of our disinfecting cleaners, do no autoclave

### Note:

- **Not intended for medical applications or patient use.**

## Variations

| Model No. | Description     | Dimensions |
|-----------|-----------------|------------|
| AB060     | Forceps, Rubber | 6"L        |
| AB061     | Forceps, Rubber | 8"L        |
| AB062     | Forceps, Rubber | 10"L       |
| AB063     | Forceps, Rubber | 12"L       |

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI