

Forceps, CutMATE

These CutMATE forceps are designed to produce consistent tissue specimen thicknesses during post-mortem preparation. The ergonomically angled, textured prongs stabilize tissue while guiding uniform 2 mm, 3 mm, or 4 mm sections sized to fit standard cassettes.



SPECIFICATIONS

Construction

- Trident-style design
- Inner prong surface: Textured
- Long ergonomically angled prongs

Performance

- Produces consistent 2 mm, 3 mm, or 4 mm tissue sections
- Helps stabilize slippery or difficult-to-hold specimens
- Fits standard plastic universal cassettes

Dimensions

- Overall length: 7"
- Prong length: 35 mm

Materials

- 304 Stainless Steel

Note:

- **Not intended for medical applications or patient use.**

Variations

Model No.	Description	Block Size
AB071	Forceps, CutMATE	Multipack, set of 4
AB072	Forceps, CutMATE	2 mm block
AB073	Forceps, CutMATE	3 mm block
AB074	Forceps, CutMATE	4 mm block
AB075	Forceps, CutMATE	2, 3, 4 mm block

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI