

Forceps, Dissection Tissue, Serrated

These serrated dissection tissue forceps are designed for versatile tissue handling during post-mortem dissection. The serrated tips provide a secure, controlled grasp for grasping, retracting, and counter-traction applications.



AB137

SPECIFICATIONS

Construction

- Tips: Serrated

Materials

- German grade stainless steel

Performance

- Supports grasping and tissue manipulation
- Suitable for retraction and counter-traction

Dimensions

- 6" (152 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI