

Forceps, Micro Grooved Tissue, Curved

These micro grooved tissue forceps are designed for precise handling and extraction during delicate post-mortem dissection. The curved, micro-fine grooved tips provide controlled access, while the unique octagonal handle offers a broad finger-grip area for precise manipulation.



AB154

SPECIFICATIONS

Construction

- Tips: Curved, micro-fine, grooved
- Handle: Octagonal with broad finger-grip area

Materials

- Premium grade stainless steel

Dimensions

- 3.75" (95 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI