

# Forceps, Extracting, #150A

These #150A extracting forceps are designed for post-mortem dental removal of upper incisors, cuspids, bicuspid, and roots. The forceps provide secure, controlled handling during extraction procedures.



AB303

## SPECIFICATIONS

### Performance

- Designed for upper incisors, cuspids, and bicuspid
- Suitable for root extraction

### Dimensions

- 7.125" (181 mm) L

### Note:

- Not intended for medical applications or patient use.

**MOPEC**  
GROUP

**Mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI