


Forceps, Potts-Smith Dressing, Tungsten Carbide

These Potts-Smith dressing forceps are designed for secure, precise tissue handling during post-mortem procedures. The serrated tips provide a controlled grasp, while the ribbed grip supports secure handling. Tungsten carbide construction provides long-lasting durability.



SPECIFICATIONS

Performance

- Designed for secure tissue handling
- Tungsten carbide construction supports long-lasting use

Construction

- Tips: Serrated
- Grip: Ribbed

Materials

- Tungsten carbide

Dimensions

- Available in 6", 7", 8.25", or 9.75" L

Note:

- **Not intended for medical applications or patient use.**

Variations

Model No.	Description	Dimensions
AB401	Forceps, Potts-Smith	6" L
AB404	Forceps, Potts-Smith	7" L
AB405	Forceps, Potts-Smith	8.25" L
AB406	Forceps, Potts-Smith	9.75" L

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI