

DISPOSABLE

Forceps, Dressing, Disposable – 12 pack

These disposable dressing forceps provide a convenient, exposure-free option for tissue and material handling during post-mortem procedures. Supplied in packs of 12, they are available in two lengths to accommodate different handling needs.



SPECIFICATIONS

Construction

- Tips: Serrated
- Grip: Ribbed

Dimensions

- Available in 4" or 4.5" L

Quantity

- 12 Pack

Materials

- Stainless steel

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Length
AB061	Forceps, Disposable	4" L
AB062	Forceps, Disposable	4.5" L

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI