

Forceps, Hemostatic, Halstead Mosquito

These Halstead mosquito hemostatic forceps are designed for secure clamping and tissue handling during post-mortem procedures. The serrated tips, finger-ring grip, and interlocking cross-arm design provide controlled clamping, with straight or curved tip options available.



SPECIFICATIONS

Construction

- Tips: Serrated
- Grip: Finger ring
- Adjustable interlocking cross arms
- Available with straight or curved tips

Materials

- Stainless steel
- Available in Standard or Premium grades

Dimensions

- 5" (127 mm) L

Note:

- **Not intended for medical applications or patient use.**

Variations

Model No.	Description	Grade	Edge
AC001	Forceps, Halstead Mosquito	Standard	Straight
AC002	Forceps, Halstead Mosquito	Standard	Curved
AC101	Forceps, Halstead Mosquito	Premium	Straight
AC102	Forceps, Halstead Mosquito	Premium	Curved

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI