

HEMOSTATIC

Forceps, Hemostatic Kelley

These Kelly hemostatic forceps are designed for secure tissue clamping during post-mortem procedures. The serrated tips, finger-ring grip, and adjustable interlocking cross arms provide firm, controlled clamping, with straight or curved tip options available.



SPECIFICATIONS

Construction

- Tips: Serrated
- Grip: Finger ring
- Adjustable interlocking cross arms
- Available with straight or curved tips

Materials

- Stainless steel
- Available in Standard or Premium grades

Dimensions

- 5.5" (140 mm) L

Note:

- **Not intended for medical applications or patient use.**

Variations

Model No.	Description	Grade	Edge
AC003	Forceps, Halstead Kelly	Standard	Straight
AC004	Forceps, Halstead Kelly	Standard	Curved
AC103	Forceps, Halstead Kelly	Premium	Straight
AC104	Forceps, Halstead Kelly	Premium	Curved

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI