

DISPOSABLE

Forceps, Hemostatic, Disposable – 12 Pack

These disposable hemostatic forceps are designed for secure clamping during laboratory and post-mortem procedures. The serrated tips provide a controlled grip, with straight or curved configurations available in two lengths.



SPECIFICATIONS

Construction

- Tips: Serrated
- Available in straight or curved configurations

Materials

- Stainless steel

Performance

- Provides secure, controlled clamping
- Suitable for a variety of laboratory and post-mortem procedures

Dimensions

- Available in 5" or 5.5" L

Quantity

- Pack of 12

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Forcep Edge	Length
AC601	Forceps, Hemostatic	Straight	5" L
AC602	Forceps, Hemostatic	Curved	5" L
AC603	Forceps, Hemostatic	Straight	5.5" L
AC604	Forceps, Hemostatic	Curved	5.5" L

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI