

Fluid Scoop, Moore Gall Stone

This Moore Gall Stone fluid scoop features an extra-long, angled handle for easy access to the abdominal and chest cavities. Its large scoop provides ample capacity for efficient fluid removal while allowing for controlled handling during autopsy procedures.



AD002

SPECIFICATIONS

Materials

- Premium grade German stainless steel

Dimensions

- 11" (279 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI