

Retractor, Blunt Flesh, 4 Prong

This curved, four prong retractor is designed for reliable retraction of exterior skin and tissue during autopsy procedures. Its specialized design provides secure tissue separation and improved access to the chest and abdominal cavities during full and limited postmortem examinations.



AD003

SPECIFICATIONS

Construction

- (4) Four prongs

Materials

- Premium grade German stainless steel

Dimensions

- 8.5" (215 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI