


Teasing Needle, Wood Handle

This teasing needle is designed for a variety of postmortem dissection and teaching applications. Its versatile design allows for precise tissue manipulation and detailed dissection procedures.



SPECIFICATIONS

Construction

- Available in a Straight or Curved point

Materials

- Wood handle with stainless steel tip

Dimensions

- AD011 & AD012: 5.5" L

- **Note:** Not intended for medical applications or patient use.

Variations

Model No.	Description	Point	Dimensions
AD011	Teasing Needle, Wood Handle	Straight	5.5" L
AD012	Teasing Needle, Wood Handle	Curved	5.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI