


# Probe, Pediatric Bowman

This Bowman probe is made of sterling silver for reliable use during postmortem dissection. Its specialized design allows for controlled probing and examination of anatomical structures.



## SPECIFICATIONS

### Construction

- Double ended

### Materials

- Sterling silver

### Dimensions

- 5" L
- AD023: 0.7 mm with 0.8 mm tips
- AD024: 0.5 mm with 0.45 mm tips

- **Note:** Not intended for medical applications or patient use.

## Variations

| Model No. | Description             | Dimensions               |
|-----------|-------------------------|--------------------------|
| AD023     | Probe, Pediatric Bowman | 0.7 mm with 0.8 mm tips  |
| AD024     | Probe, Pediatric Bowman | 0.5 mm with 0.45 mm tips |

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI