

Retractor/Lifter, Sharp/Blunt, Curved

This curved sharp/blunt retractor and lifter features a paddle end designed for isolating and retracting tissue during postmortem dissection. Its specialized design provides controlled tissue manipulation and improved access to anatomical structures during detailed autopsy procedures.



SPECIFICATIONS

Construction

- Three prongs
- Available in Sharp or Blunt prong type

Materials

- Stainless steel

Dimensions

- AD030 & AD031: 6.5" (165 mm) L

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Prong Type	Dimensions
AD030	Retractor/Lifter, Curved	Sharp	6.5" L
AD031	Retractor/Lifter, Curved	Blunt	6.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI