

Retractor, Blunt/Sharp

This high-quality retractor is designed to separate and hold tissue during autopsy procedures for improved access and visibility. Its sturdy design allows for deeper, more controlled cutting and probing during detailed postmortem dissection.



SPECIFICATIONS

Construction

- Available in Sharp or Blunt prong type
- 3 X 4 prongs

Materials

- Stainless steel

Dimensions

- AD047 & AD048: 5.5" L
- AD049 & AD050: 9.5" L
- 5.5" L has a 2.75" spread
- 9.5" L has a 4.75" spread

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Prong Type	Dimensions	Spread
AD047	Retractor, Blunt/Sharp	Sharp	5.5" L	2.75"
AD048	Retractor, Blunt/Sharp	Blunt	5.5" L	2.75"
AD049	Retractor, Blunt/Sharp	Sharp	9.5" L	4.75"
AD050	Retractor, Blunt/Sharp	Blunt	9.5" L	4.75"

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI